


Software Requirements

Pre-operative Planning Software

Your planning helper for surgeries



Service

- Standalone software installed on a single PC;
- Internet connection is required;
- Periodic update of implant database.

Minimum system requirements

- Processor: x86 and compatible with CPU (Intel/AMD), min. 2.0 GHz, dual core;
- RAM: min. 4GB;
- Graphics card - NVIDIA GeForce 9xxx, 650 GT, ATI Radeon HD R5 Series, with support OpenGL 4.0;
- Graphics memory: min. 512 MB;
- Disk space: 5 GB of free disk space;
- MacOSX® High Sierra or later version;
- Microsoft Windows® 7 or later version.

Recommended system requirements

- Processor: x64 and compatible with CPU (Intel/AMD), min. 2.5 GHz, i5 or higher;
- RAM: 8GB or higher;
- Graphics card: OpenGL 4.0 support;
- Graphics memory: min. 1GB DDR5;
- Disk space: 5 GB of free disk space;
- MacOSX® High Sierra or later version;
- Microsoft Windows® 7 or later version;
- An antivirus to help protecting the images retrieved from the PACS.

Contact us

For more information or any question please visit www.peekmed.com or send an email to info@peekmed.com.

